

ATCC Medium: 1827 BHI w/ heat-inactivated Fetal Bovine Serum and Glucose

Base Medium

Brain Heart Infusion (BD 237500).....37.0 g
DI Water.....850.0 ml

Autoclave at 121°C.

Cool to 50°C and aseptically add the following:

Additive

Heat inactivated fetal bovine serum.....100.0 ml
Glucose Solution (see below).....50 ml

Glucose Solution

Glucose.....2.0 g
DI Water.....50 ml

Filter sterilize.

Prepare medium anaerobically using 80% N₂, 10% CO₂, 10% H₂.