

ATCC Medium: 234 Peptone-Succinate Agar

(NH ₄) ₂ SO ₄	1.0 g
MgSO ₄ x 7H ₂ O.....	1.0 g
MnSO ₄ x H ₂ O.....	2.0 mg
FeCl ₃ x 6H ₂ O.....	2.0 mg
Succinic Acid.....	1.68 g
Peptone.....	5.0 g
Agar.....	1.5 g
DI Water.....	1000 ml

Adjust pH to 7.0. Autoclave at 121°C.