

ATCC medium: 2261 Modified R2A medium

Yeast extract.....	0.5 g
Proteose Peptone No. 3 (BD 211693).....	0.5 g
Casamino acids.....	0.5 g
Glucose.....	0.5 g
Tween 80	0.1% (v/v)
Sodium pyruvate.....	0.3 g
K ₂ HPO ₄	0.3 g
MgSO ₄ . 7H ₂ O	0.05 g
Agar.....	15.0 g
Distilled deionized water.....	1.0 L

Adjust pH to 7.2 with potassium phosphates before adding agar.
Autoclave at 121C for 15 minutes.