

**ATCC medium: 2084 *Halomonas* medium**

NaCl.....	50.0 g
Casamino Acids (with Vitamins).....	7.5 g
Proteose Peptone No. 3 (BD 211693).....	5.0 g
Yeast extract.....	1.0 g
Trisodium citrate.....	3.0 g
MgSO <sub>4</sub> . 7H <sub>2</sub> O .....	20.0 g
K <sub>2</sub> HPO <sub>4</sub> .....	0.5 g
Fe(NH <sub>4</sub> ) <sub>2</sub> (SO <sub>4</sub> ) <sub>2</sub> . 6H <sub>2</sub> O .....	0.05 g
Agar (if necessary).....	15.0 g
Distilled water.....	1.0 L

Adjust pH to 7.3 with KOH.

Autoclave at 121C for 15 minutes.