

ATCC medium: 1233 Revised anthranilic acid medium for ATCC 31622

Na₂HPO₄6.0 g
KH₂PO₄3.0 g
NaCl.....0.5 g
NH₄Cl1.0 g
Distilled water.....950.0 ml

Autoclave basal medium for 15 minutes at 121C. Cool and aseptically add the following sterile solutions:

0.01 M CaCl₂10.0 ml
0.1 M MgSO₄10.0 ml
20% Glucose.....25.0 ml
1% Anthranilic acid (in ethanol)....5.0 ml