

**ATCC medium: 2115 *Halanaerobium lacusrosei* medium**

Yeast extract.....10.0 g  
Peptone.....2.0 g  
NaCl.....200.0 g  
MgSO<sub>4</sub>.....0.5 g  
K<sub>2</sub>HPO<sub>4</sub>.....1.0 g  
Agar.....15.0 g  
Distilled water.....1.0 L

Adjust pH to 7.1. Boil under N<sub>2</sub> gas. Dispense under gas.  
Autoclave for 15 minutes at 121C.