

**ATCC MEDIUM: 1187 Mannan Salts Medium**

Yeast Mannan	3.0	g
(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	0.5	g
MgSO <sub>4</sub>	0.2	g
FeSO <sub>4</sub>	0.01	g
CaCl <sub>2</sub> x H <sub>2</sub> O	0.05	g
Yeast Extract	0.5	g
K <sub>2</sub> HPO <sub>4</sub>	7.54	g
KH <sub>2</sub> PO <sub>4</sub>	2.32	g
DI Water	1000	ml

*Autoclave at 121 C for 15 minutes.*