

ATCC medium: 1273 LB medium (ATCC medium 1065) with 20 mcg/ml tetracycline

LB medium (ATCC Medium 1065, see below) with 20.0 mcg/ml tetracycline

ATCC Medium 1065:

| | |
|--------------------------------|--------|
| Tryptone (BD 211705)..... | 10.0 g |
| Yeast Extract (BD 212750)..... | 5.0 g |
| NaCl..... | 10.0 g |
| Distilled water..... | 1.0 L |

Autoclave at 121C for 15 minutes.