

**ATCC medium: 1210 Avian mycoplasma medium**

PPL0 Broth w/o CV (BD 255420).....	14.7 g
Distilled or deionized water.....	650.0 ml
Swine or horse serum (heat-inactivated).....	150.0 ml
Yeast Extract Solution (GIBCO 18180).....	100.0 ml
Glucose solution [10% (wt/vol)].....	10.0 ml
Arginine solution [10% (wt/vol)].....	10.0 ml
NAD [1% (wt/vol)].....	10.0 ml
Phenol red solution [0.1% (wt/vol)] (for liquid medium only)....	20.0 ml
Agar (for solid medium only)*.....	10.0 g

Mix the PPL0 broth base in the water (together with agar for solid medium) and autoclave at 121C for 15 minutes. Mix the remaining ingredients and heat to 56C in a water bath. Combine both lots of ingredients, mix, and dispense.

\*Use agar known to be noninhibiting to mycoplasmas.