

ATCC Medium 4: Rabbit Blood Medium

Agar Medium

Blood Agar Base (BD 211037)..... 40 g
DI Water..... 950 mL

Broth Medium

Heart Infusion Broth..... 2 g
Pancreatic Digest of Casein..... 13 g
Yeast Extract..... 5 g
Sodium Chloride (NaCl)..... 5 g
DI Water 950 mL

Adjust pH to 7.3 ± 0.2 (broth)

Mix medium thoroughly and autoclave at 121° C.

Cool to 48-50°C and aseptically add the following:

Rabbit Blood..... 50 mL

Gently mix and dispense into appropriate vessel.