

**ATCC Medium 0024: Nutrient Agar/Broth w/ 1% Potato Starch**

**Agar Base Medium**

Nutrient Agar (BD 213000).....	23.0 g
Potato Starch.....	10.0 g
DI Water.....	1000.0 mL

**Broth Base Medium**

Nutrient Broth (BD 234000).....	8.0 g
Potato Starch.....	10.0 g
DI Water.....	1000.0 mL

Adjust pH to **6.8 +/- 0.2**. Autoclave at 121°C.

Attempt to keep the starch in suspension when dispensing.