

ATCC medium: 1760 *Veillonella* medium (DSM)

Trypticase Peptone (BD 211921).....5.0 g
Yeast extract.....3.0 g
Sodium (DL) lactate.....7.5 g
Sodium thioglycollate.....0.75 g
Tween 80.....1.0 g
Glucose.....1.0 g
Putrescine.....3.0 mg
Resazurin.....1.0 mg
Distilled water.....1.0 L

Adjust medium to pH 7.5 with solid K_2CO_3 . Prepare and dispense medium anaerobically under 100% N_2 atmosphere. Autoclave at 121C for 15 minutes.