



## Product Sheet

# pAB366Ben (ATCC® 77136™)

Please read this FIRST



Storage Temp.  
**>Store unopened  
frozen vial at -  
80°C.**

---



Biosafety Level  
**1**

### Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

### Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: pAB366Ben (ATCC® 77136™)

### Shipping Information

Frozen glycerol stock of *Escherichia coli* containing the plasmid.

American Type Culture Collection  
PO Box 1549  
Manassas, VA 20108 USA  
[www.atcc.org](http://www.atcc.org)

800.638.6597 or 703.365.2700  
Fax: 703.365.2750  
Email: [Tech@atcc.org](mailto:Tech@atcc.org)

Or contact your local distributor

## Description

This is a shuttle reporter plasmid for use in *Aspergillus nidulans* or *Neurospora crassa* that can be used to test expression of bleomycin resistance in transformants. Expresses the mutant beta-tubulin gene conferring benomyl resistance. The order of the major features in this plasmid is : pUC8 EcoRI BamHI am promoter ble XhoI am 3 region tub-2 (benomyl resistance) pUC8.A 2.5 kb Sall/HindIII fragment contains the complete tubulin gene and 5 and 3 flanking regions. Restriction digests of the vector gave the following sizes (kb): BamHI 5.1, 2.1; HindIII 4.0, 3.2; BamHI/HindIII 2.6, 2.5, 1.5, 0.65; BamHI/XhoI 3.0, 1.35, 1.2, 0.9, 0.8; Sall/HindIII 3.3, 2.6, 0.85, 0.55. -----ATCC staff

**Designation:** pAB366Ben

## Propagation

Transfer a loopful to a test tube containing 5 mL LB+50mg/mL of ampicillin broth. A loopful of culture can also be streaked on an LB + amp agar plate. Incubate cultures at 37°C. Isolate DNA using standard plasmid preparation procedures.

### Growth Conditions

**Temperature:** 37.0°C

### Medium

ATCC® Medium 1227: LB Medium (ATCC medium 1065) with 50 mcg/ml ampicillin

## Vector Information

Size (kb): 7.1999998092651370

Vector: pAB366Ben (plasmid)

Construction: pAB366, pBT6

Marker(s): ampR, bleR, benR

Construct size (kb): 7.199999809265137

Features: marker(s): ampR, bleR, benR

replicon: pMB1

## References

References and other information relating to this product are available online at [www.atcc.org](http://www.atcc.org).

## Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the current publication of the *Biosafety in Microbiological and Biomedical Laboratories* from the U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes for Health.

## ATCC Warranty

The viability of ATCC® products is warranted for 30 days from the date of shipment, and is valid only if the product is stored and cultured according to the information included on this product information sheet. ATCC lists the media formulation that has been found to be effective for this strain. While other, unspecified media may also produce satisfactory results, a change in media or the absence of an additive from the ATCC recommended media may affect recovery, growth and/or function of this strain. If an alternative medium formulation is used, the ATCC warranty for viability is no longer valid.

## Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. ATCC is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of cultures.



Product Sheet

## pAB366Ben (ATCC® 77136™)

Please see the enclosed Material Transfer Agreement (MTA) for further details regarding the use of this product. The MTA is also available on our Web site at [www.atcc.org](http://www.atcc.org)

Additional information on this culture is available on the ATCC web site at [www.atcc.org](http://www.atcc.org).  
© ATCC 2013. All rights reserved. ATCC is a registered trademark of the American Type Culture Collection. [02/06]

Please read this FIRST

The graphic is a black rectangle with white text and icons. At the top left is a thermometer icon. To its right, the text reads: "Storage Temp. >Store unopened frozen vial at -80°C." Below this is a horizontal dotted line. At the bottom left is a biohazard symbol. To its right, the text reads: "Biosafety Level 1".

### Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

### Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: pAB366Ben (ATCC® 77136™)

### Shipping Information

Frozen glycerol stock of *Escherichia coli* containing the plasmid.

American Type Culture Collection  
PO Box 1549  
Manassas, VA 20108 USA  
[www.atcc.org](http://www.atcc.org)

800.638.6597 or 703.365.2700  
Fax: 703.365.2750  
Email: [Tech@atcc.org](mailto:Tech@atcc.org)

Or contact your local distributor