




## Product Sheet


# pNKY70 (ATCC® 77183™)

Please read this **FIRST**



Storage Temp.  
**Store unopened  
freeze dried vial  
at 4°C.**

---



Biosafety Level  
**1**

### Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

### Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: pNKY70 (ATCC® 77183™)

### Shipping Information

Freeze dried *E. coli* containing the plasmid.

American Type Culture Collection  
PO Box 1549  
Manassas, VA 20108 USA  
[www.atcc.org](http://www.atcc.org)

800.638.6597 or 703.365.2700  
Fax: 703.365.2750  
Email: [Tech@atcc.org](mailto:Tech@atcc.org)

Or contact your local distributor

## Description

This is a YI-type (yeast integrating) vector that can be used for promoter/enhancer trap screens. It contains a HIS4C/URA3 fusion cassette that can be used to make fusions to other proteins of interest. It also includes an M13 sequencing primer approximately 300 bp from the HIS4C/URA3 fusion junction. It encodes 93 amino acids of HIS4C and the complete 267 amino acids of URA3. This vector was constructed from pNKY48 (ATCC#77182) by inserting an EcoRI 8-mer at the HindIII site in the polylinker. The order of the major features in this plasmid is: M13 primer ClaI XhoI Sall XbaI (stop) BamHI EcoRI HIS4C bridge/URA3 BamHI HindIII Bgl II NcoI EcoRI M13 ori pMB1 ori ampR. - Genetics 117: 5-12, 1987

**Designation:** pNKY70

## Propagation

1. Open vial according to instructions. 2. Aseptically add 0.3 to 0.4 mL of liquid medium to the freeze-dried pellet and mix well. Transfer 100  $\mu$ L to a test tube containing 5 mL LB+ ampicillin (50-100  $\mu$ g/mL). A loopful of culture can also be streaked on an agar plate of the same. Incubate cultures at 37°C. 3. Isolate DNA using standard plasmid preparation procedures.

### Growth Conditions

**Temperature:** 37.0°C

### Medium

ATCC® Medium 1227: LB Medium (ATCC medium 1065) with 50 mcg/ml ampicillin

## Vector Information

Size (kb): 4.4000000953674320

Vector: pNKY70 (plasmid)

Construction: pNKY48 (ATCC 77182)

Marker(s): URA3, ampR

Construct size (kb): 4.400000095367432

Features: marker(s): ampR, URA3

replicon: pMB1, M13

terminator: URA3

## References

References and other information relating to this product are available online at [www.atcc.org](http://www.atcc.org).

## Notes

Restriction digests of the vector gave the following sizes (in kb): Bgl II 4.4; EcoRI 2.9, 1.5; PstI 4.4; BamHI 3.2, 1.2; XbaI 4.4. ATCC Staff

## Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the current publication of the *Biosafety in Microbiological and Biomedical Laboratories* from the U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes for Health.

## ATCC Warranty

The viability of ATCC® products is warranted for 30 days from the date of shipment, and is valid only if the product is stored and cultured according to the information included on this product information sheet. ATCC lists the media formulation that has been found to be effective for this strain. While other, unspecified media may also produce satisfactory results, a change in media or the absence of an additive from the ATCC recommended media may affect recovery, growth and/or function of this strain. If an alternative medium formulation is used, the ATCC warranty for viability is no longer valid.

## Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product




Product Sheet

## pNKY70 (ATCC® 77183™)


---

### Please read this FIRST



Storage Temp.  
**Store unopened  
freeze dried vial  
at 4°C.**

---



Biosafety Level  
**1**

---

### Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

---

### Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: pNKY70 (ATCC® 77183™)

---

### Shipping Information

Freeze dried *E. coli* containing the plasmid.

sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. ATCC is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of cultures.

Please see the enclosed Material Transfer Agreement (MTA) for further details regarding the use of this product. The MTA is also available on our Web site at [www.atcc.org](http://www.atcc.org)

Additional information on this culture is available on the ATCC web site at [www.atcc.org](http://www.atcc.org).

© ATCC 2013. All rights reserved. ATCC is a registered trademark of the American Type Culture Collection. [02/06]

American Type Culture Collection  
PO Box 1549  
Manassas, VA 20108 USA  
[www.atcc.org](http://www.atcc.org)

800.638.6597 or 703.365.2700  
Fax: 703.365.2750  
Email: [Tech@atcc.org](mailto:Tech@atcc.org)

Or contact your local distributor