





Product Sheet

p413 ADH (ATCC® 87370™)

Please read this FIRST



Storage Temp.
**Store unopened
freeze dried vial
at 4° C.**



Biosafety Level
1

Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: p413 ADH (ATCC® 87370™)

Shipping Information

Freeze dried *E. coli* containing the plasmid

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor

Description

This is a low copy number shuttle expression vector. One of 32 yeast expression vectors (ATCC#87354-87385) differing in promoter (GPD, TEF, ADH, CYC1), selection marker (HIS3, TRP1, LEU2, URA3) and replicon (CEN6/ARSH4, 2 micron). This plasmid, p413 ADH, has the ADH promoter, the HIS3 marker and the CEN6/ARSH4 replicon. The wild type ADH promoter is active when cells are grown in glucose media but can be repressed 2-10 fold on non-fermentable carbon sources.

Designation: p413 ADH

Propagation

1. Open vial according to instructions. 2. aseptically add 0.3 to 0.4 mL of liquid medium to the freeze-dried pellet and mix well. Transfer 100 μ L to a test tube containing 5 mL LB+ ampicillin (50-100 μ g/mL). A loopful of culture can also be streaked on an agar plate of the same. Incubate cultures at 37°C. 3. Isolate DNA using standard plasmid preparation procedures.

Growth Conditions

Temperature: 37.0°C

Medium

ATCC® Medium 1227: LB Medium (ATCC medium 1065) with 50 mcg/ml ampicillin



Vector Information

Size (kb): 6.696

Vector: p413 ADH (plasmid)

Promoters: Promoter for expression ADH

Construction: pRS413

Marker(s): ampR, HIS3

Construct size (kb): 6.696

Features: marker(s): HIS3

marker(s): ampR

promoter for expression: ADH

replicon: CEN6/ARSH4

MCS: XhoI...XbaI

terminator: CYC1



References

References and other information relating to this product are available online at www.atcc.org.



Notes

Restriction digests of the clone gave the following sizes (in kb): BamHI 6.8 ; XhoI 6.8 ;
SacI/XbaI 5.2, 1.6. ATCC Staff



Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the current publication of the *Biosafety in Microbiological and Biomedical Laboratories* from the U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes for Health.

ATCC Warranty


The viability of ATCC® products is warranted for 30 days from the date of shipment, and is valid only if the product is stored and cultured according to the information included on this product information sheet. ATCC lists the media formulation that has been found to be effective for this strain. While other, unspecified media may also produce satisfactory results, a change in media or the absence of an additive from the ATCC recommended media may affect recovery, growth and/or function of this strain. If an alternative medium formulation is used, the ATCC warranty for viability is no longer valid.




Product Sheet

p413 ADH (ATCC® 87370™)

Please read this FIRST



Storage Temp.
**Store unopened
freeze dried vial
at 4⁰ C.**



Biosafety Level
1

Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: p413 ADH (ATCC® 87370™)

Shipping Information

Freeze dried *E. coli* containing the plasmid

Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. ATCC is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of cultures.

Please see the enclosed Material Transfer Agreement (MTA) for further details regarding the use of this product. The MTA is also available on our Web site at www.atcc.org

Additional information on this culture is available on the ATCC web site at www.atcc.org.

© ATCC 2017. All rights reserved. ATCC is a registered trademark of the American Type Culture Collection. [12/08]

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor