





Product Sheet

pCM190 (ATCC® 87662™)

Please read this FIRST



Storage Temp.
**Store unopened
frozen vial at -80⁰
C.**



Biosafety Level
1

Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: pCM190 (ATCC® 87662™)

Shipping Information

Frozen glycerol stock of *Escherichia coli* DH5alpha containing the plasmid.

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor

Description

This is a yeast episomal (YE-type) expression vector driven by the tetO-CYC1 promoter. This vector is one of 7 tetracycline-regulated expression vectors (ATCC#87656 87662) differing in the level of expression, mode of replication and selectable marker. Restriction digests of the vector gave the following sizes (in kb): BamHI 8.0; PstI 8.0; BamHI/XhoI 6.9, 0.8. ---

-ATCC staff

Designation: pCM190

Propagation

Transfer a loopful to a test tube containing 5 mL LB+50mg/mL of ampicillin broth. A loopful of culture can also be streaked on an LB + amp agar plate. Incubate cultures at 37⁰C. Isolate DNA using standard plasmid preparation procedures.

Growth Conditions

Temperature: 37.0°C

Medium

ATCC® Medium 1227: LB Medium (ATCC medium 1065) with 50 mcg/ml ampicillin

Vector Information

Size (kb): 8.0120000839233400

Vector: pCM190 (plasmid)

Promoters: Promoter CMV

Marker(s): URA3, ampR

Construct size (kb): 8.01200008392334

Features: marker(s): URA3

marker(s): ampR

operator: tetO - seven repeats

other: tTA transcriptional activator

promoter: CMV

promoter: CYC1

replicon: 2 micron

replicon: pMB1

MCS: BamHI...PstI

terminator: ADH1

terminator: CYC1

References

References and other information relating to this product are available online at www.atcc.org.

Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the current publication of the *Biosafety in Microbiological and Biomedical Laboratories* from the U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes for Health.

ATCC Warranty

The viability of ATCC® products is warranted for 30 days from the date of shipment, and is valid only if the product is stored and cultured according to the information included on this product information sheet. ATCC lists the media formulation that has been found to be effective for this strain. While other, unspecified media may also produce satisfactory results, a change in media or the absence of an additive from the ATCC recommended media may affect recovery, growth and/or function of this strain. If an alternative medium formulation is used, the ATCC warranty for viability is no longer valid.

Disclaimers


This product is intended for laboratory research purposes only. It is not intended for use in




Product Sheet

pCM190 (ATCC® 87662™)

Please read this FIRST



Storage Temp.
**Store unopened
frozen vial at -80°
C.**



Biosafety Level
1

Intended Use

This product is intended for research use only. It is not intended for any animal or human therapeutic or diagnostic use.

Citation of Strain

If use of this culture results in a scientific publication, it should be cited in that manuscript in the following manner: pCM190 (ATCC® 87662™)

Shipping Information

Frozen glycerol stock of *Escherichia coli* DH5alpha containing the plasmid.

humans.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. ATCC is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of cultures.

Please see the enclosed Material Transfer Agreement (MTA) for further details regarding the use of this product. The MTA is also available on our Web site at www.atcc.org

Additional information on this culture is available on the ATCC web site at www.atcc.org.

© ATCC 2013. All rights reserved. ATCC is a registered trademark of the American Type Culture Collection. [02/06]

American Type Culture Collection
PO Box 1549
Manassas, VA 20108 USA
www.atcc.org

800.638.6597 or 703.365.2700
Fax: 703.365.2750
Email: Tech@atcc.org

Or contact your local distributor